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2024 Zhonglun Junior High School Science Fair Challenge Activity

Checkpoint Name: Determination of Reducing Sugars

Objective of the Challenge: Students will deepen their understanding of the basic concepts of the junior high school biology nutrition unit through hands-on experiments. They will also attempt to explain and apply these concepts, fostering a spirit of scientific inquiry.

Challenge Process: Conduct the experiment and fill in the results on the worksheet. The challenge is completed once all questions are answered.

Checkpoint Content:

Checkpoint: Add 2.5ML of each of the following liquids into separate dropper bottles: Super Supau sports drink, Coca-Cola, Coca-Cola Zero, and Pomelo vegetable juice. Then, add 2.5ML of Benedict's solution to each. Determine which liquid is the sweetest based on the final color. (Sweetness increases from sugar-free to the sweetest in the following order: blue → green → yellow → orange → red).

Materials Needed for the Checkpoints: :

1. Test tubes X20
2. Test tube racks X5
3. Cloths X3

- 4.3ML pipettes X5
- 5.10ML graduated cylinders X5
- 6.5ML dropper bottles X400
- 7.250ML beakers X5
- 8.1000W induction cooker X1
- 9.Supau sports drink 1L X2
- 10.Coca-Cola 2L X1
- 11.Coca-Cola Zero 2L X1
12. Bomi vegetable juice 250ML X8
- 13.Benedict' s solution 500ML X4
- 14.Inner pot for electric cooker 4L X1
- 15.Fine-tip oil-based pens X5
- 16.Label paper
- 17.Approval stamp 5% starch solution 500ML

4. Core Curriculum Content Corresponding to the Checkpoints:

- **PE-IV-2:**Ability to correctly and safely operate appropriate items, equipment, instruments, technological devices, and resources for the learning stage. Ability to conduct objective qualitative observations or quantitative measurements and record them accurately.
- **AI-IV-3:**Using learned scientific knowledge and various methods of scientific exploration to explain the causes of natural phenomena, building confidence in scientific learning.

5.Script Process:

1. Pre-experiment Testing and Training:

- Conduct pre-experiment tests to determine the final color results of four types of beverages.
- Train the reception students (divided into 3 groups: Group 1 3persons and Group 2 3persons each get two test tube racks, Group 3 (bilingual group) gets one test tube rack).

2. Preparation of Dropper Bottles:

- On the day before the science fair, pre-fill the dropper bottles: label the bottles with numbers 1, 2, 3, and 4 using a marker.
- Fill bottle 1 with 2.5ML of Supau sports drink, bottle 2 with 2.5ML of Coca-Cola, bottle 3 with 2.5ML of Coca-Cola Zero, and bottle 4 with 2.5ML of Pomi mixed vegetable and fruit juice.

3. Preparation of Test Tube Racks:

- Prepare 5 test tube racks; place 4 test tubes in each rack (label the test tubes with numbers 1, 2, 3, and 4 corresponding to the dropper bottles).

4. Distribution to Participants:

- Reception students distribute one set of test tubes and a test tube rack to each participant, along with dropper bottles numbered 1, 2, 3, and 4, and a 3ML pipette.

5. Guidance for Participants:

- Reception students guide participants to squeeze the liquid from the dropper bottles into the corresponding test tubes (label the test tubes with numbers).
- After this, use the pipette to draw 2.5ML of Benedict' s solution into the first test tube. Remind participants to check the pipette for the 0.5ML to 3.0ML scale.
- After drawing 2.5ML of liquid into the pipette, squeeze it into a graduated cylinder to check if it is also 2.5ML. If it matches, add it directly to the corresponding test tube. For the third and fourth test tubes, directly draw 2.5ML of liquid with the pipette, skipping the graduated cylinder measurement step.

6. Water Bath Heating:

- Participants place the four test tubes into a beaker containing 100ML of hot water (be careful when getting hot water from the dispenser!).
- Place the beaker into a 500ML water bath for heating, using an induction cooker set to medium power (around 400W) to maintain the temperature. Ensure the beaker does not float and tip over.
- Participants who complete the above steps can place their beakers into the water bath, with up to five beakers at a time.

7. Observation and Worksheet Completion:

- Participants observe the results (color) of the four test tubes and complete the worksheet, comparing which beverage is the sweetest.

8. Stamping the Worksheet:

- Reception students stamp the worksheets of participants who have completed them.

2024 Zhonglun Junior High School Science Fair Challenge Activity

Checkpoint Worksheet

Subject: Determination of Reducing Sugars

School: _____ Grade: _____ Gender: _____ Name: _____


Sugar is not a necessity for the body and can even be considered harmful. Its negative effects include: increasing body fat, causing metabolic syndrome, dulling brain function, hindering work or daily life, inducing anxiety, leading to an unstable mental state, and triggering various diseases such as diabetes, hypertension, arteriosclerosis, allergies, increasing cancer risk, accelerating aging, and increasing the likelihood of dementia.

There are many types of beverages in daily life, and those you think have low sugar content often contain more sugar than expected. The main goal of this activity is to recognize hidden high-sugar foods, learn to read food nutrition labels, identify hidden sugars, and consciously control and choose your intake.

In this activity, the color change of Benedict's solution will be used to show you that beverages are actually full of hidden dangers. Even drinks that are considered healthy might be better consumed in moderation. The original color of Benedict's solution is light blue. After reacting with heated beverages, the sweetness concentration increases from sugar-free to high sugar, with color changes from blue → green → yellow → orange → red (red indicates extremely high sweetness).

Principle of Benedict's Solution Reaction:

Benedict's solution reacts with aldehydes to produce a red precipitate of cuprous oxide (Cu₂O). Glucose and galactose are aldoses, while fructose is a ketose, but its adjacent C atom has an OH- group that can convert into an aldehyde group, thus it also reacts.

Reducing sugar → Organic acid; Cu²⁺ → Cu⁺ (red cuprous oxide Cu₂O)

Oxidation

Reduction

Checkpoint Questions: (Please check the result)

Final color Supau Drink : Blue Green Yellow Orange Red

Final color Coca-Cola : Blue Green Yellow Orange Red

Final color Coke Zero : Blue Green Yellow Orange Red

Final color Bomi Juice : Blue Green Yellow Orange Red

Which drink is the sweetest? Supau Sports Drink Coca-Cola

Coke Zero Bomi Vegetable Juice

Which drink has the lowest sugar content? Supau Sports Drink

Coca-Cola Coke Zero Bomi Vegetable Juice

Write down your findings and thoughts: _____